

# Hiperskin<sup>®</sup> Fixing Rolls & Dressings

## TYPE 1.

### Hiperskin<sup>®</sup> TR!



Translucent film roll-type bandage with excellent adhesive and water-proof qualities

- Translucent material provides easy visualization of dressing
- Blocks water and germs but is permeable to air, to maintain appropriate humidity
- Excellent adhesiveness prevents falling off during showers

## TYPE 2.

### Hiperskin<sup>®</sup> NV!



Non-woven roll-type bandage with fine ventilation pores that facilitate skin respiration

- Non-woven fabric provides excellent elasticity and ventilation, allowing ventilation even after long hours of use
- Prevents sweat accumulation on skin, and is appropriate for use on areas with high-activity levels




# Hiperskin<sup>®</sup> Fixing Rolls & Dressings

## Hiperskin Premium & Hiperskin Premium Plus

### “Fix-roll dressing with excellent elasticity”

- Polyurethane nonwoven fabric of melt-blown method is used, and it is a premium fix roll with excellent elasticity, waterproofness and air permeability
- It can be used conveniently in areas where it is difficult to apply dressing such as bending area and joints.
- You can cut and use as much as you need.
- It consists of the following two types of products depending on the adhesive property.

### TYPE 3. Hiperskin<sup>®</sup> Premium!



### TYPE 4. Hiperskin<sup>®</sup> Premium<sup>Plus</sup>!



#### Hiperskin<sup>®</sup> Premium!

- It is a fix roll dressing using acrylic adhesives with excellent adhesiveness.
- It has good adhesion moment and sticks to any part.
- It does not fall well with frequent movements.

#### Hiperskin<sup>®</sup> Premium<sup>Plus</sup>!

- Hypoallergenic fix roll dressing with soft silicone adhesive.
- Excellent durability of adhesive, so adhesive strength is maintained for a long time.
- It is easy to use because of excellent handling when using the product.
- When removing the product, it is smoothly removed without the pain of the attachment surface.