

Hiperskin[®]NV

Product Name	Product Code	Size(cm x M)	Size (Roll / 1Box)
Hiperskin NV	NV2510ZA	2.5cm x 10M	1 Roll / 1 Box
	NV0510ZA	5cm x 10M	1 Roll / 1 Box
	NV1010ZA	10cm x 10M	1 Roll / 1 Box
	NV1510ZA	15cm x 10M	1 Roll / 1 Box



- Non-woven fabric roll bandage with fine ventilation pores that facilitate skin respiration
- Manufactured using non-woven fabric, which is soft and has excellent surface texture
- Excellent adhesiveness and elasticity prevent falling off even when applied to contoured areas

Hiperskin[®]TR

Product Name	Product Code	Size (cm x M)	Packing (Roll / 1Box)
Hiperskin TR	TR0510ZA	5cm x 10M	1 Roll / 1 Box
	TR1010ZA	10cm x 10M	1 Roll / 1 Box



- Translucent polyurethane film dressing material allows easy observation of inside
- Blocks water and bacteria while being permeable to air, and maintains appropriate humidity
- Skin-like softness and has excellent surface texture
- Is not easily distinguishable when applied to skin, and is aesthetically favorable
- Excellent adherence and is appropriate for showers and during swimming sessions

Hiperskin[®]Premium

Product Name	Product Code	Size (cm x M)	Packing(Roll / 1Box)
Hiperskin Premium	HR2503ZA	2.5cm x 3M	1 Roll / 1 Box
	HR0503ZA	5cm x 3M	1 Roll / 1 Box
	HR1003ZA	10cm x 3M	1 Roll / 1 Box
	HR2510ZA	2.5cm x 10M	1 Roll / 1 Box
	HR0510ZA	5cm x 10M	1 Roll / 1 Box
	HR1010ZA	10cm x 10M	1 Roll / 1 Box



- Fixing roll-type product which is secondary dressing material. For use in fixing after applying appropriate dressing material such as foam dressing material
- Does not fall off when applied to contoured surfaces or joints with wide ranges of movement, due to excellent elasticity and adhesiveness
- Soft touch and texture comforts skin with airtight sealing
- Does not generate fluff because it is made of non-woven fabric, and maintains cleanness after dressing.
- Is not easily visible when applied on skin, and is aesthetically favorable

Hiperskin[®] NV, TR, Premium

Roll-type bandage that wraps comfortable around skin with excellent elasticity

■ Appearance

Medical support aid for protection, adhesion, along with fixing of roll-type bandage, and dressing (NV, TR)

Milk-white or yellowish polyurethane film bandage with adhesive applied on one side (Premium)

■ Efficacy & Effectiveness

Prevention of infection and protection of wound site

Medical support aid for protection, adhesion, along with fixing of roll-type bandage, and dressing

■ How to Use

- 1) Wash gently and thoroughly with flowing water or saline
- 2) Dry wound site
- 3) Cut product to fit wound-site
- 4) Remove plaid peel and apply tightly on wound

■ Cautions

- 1) Check for debris on site of application. Stop using immediately in case of skin hypersensitivity
- 2) Store in dark, dry and cool area

■ Storage methods and expiration

Store at room temperature (1-30°C), airtight container, expires 36months from date of manufacture