


Hydrocolloid Dressing

Hydrocolloid dressing

“Hydrocolloid dressing to help protect and heal wounds”

Renoderm[®]




Spot type

- The hydrocolloid component absorbs the exudate to form a gel and safely decompose the necrotic tissue.
- The formed gel keeps the wound surface as moist environment and helps wound healing.
- It has good adhesion and it can be attached to the curved area. No secondary dressing is required.
- The gel formed by absorbing the exudate from the wound surface makes hydrocolloid lose adhesion and minimizes the damage of the wound when removing the product

Functional hydrocolloid dressing

ILDONG Meditouch^H Sun Protection



<자외선차단시험성적서>

“UV protection hydrocolloid”

- UV Protection Factor : UPF 50+ (More than 98% UV Protection)
- Prevents secondary damage from UV rays from outside activity.



SOLSOLPLUS



“Containing Titrated Extract of Centella Asiatica”

- Centella extract helps to moisturize and heal wounds.



56, 1gongdan-ro 5-gil, Gumi-si,
Gyeongsangbuk-do, Republic of Korea

